

Recovery

Balance

Mobility

Health



AlterG – The New Standard of Care in Rehab and Conditioning

Get patients better faster: Pain free movement encourages earlier walking with normal gait mechanics. Precise partial weight bearing for earlier incorporation of closed-kinetic chain activities.

Get patients better safely: Patients are comfortably secure inside the AlterG. Gentle air support allows for increased confidence and safe progression of activities to improve neuro-muscular control and balance.

Decrease overall cost of care: Easy to use and maintain, with minimal resource requirements from facility and staff.

Create facility revenue: Mobilizing patients faster can lead to a placement in a higher Resource Utilization Group (RUG). Attract more doctor referrals and patients with state-of-the-art NASA technology.

Be competitive: Set yourself apart and provide the best technology at your facilities.





"We love our AlterG. It serves our population perfectly. It gives us the opportunity to have our clients walk distances and at a pace they wouldn't dream of. It doesn't cure the common cold, but it sure cures the common knee or hip pains!"

Maryann Bowles, PT
 Rehab Services Manager
 Life Care Center of Littleton, CO

FDA-cleared and recommended by top physical therapists and doctors

The FDA-cleared and patented AlterG Anti-Gravity Treadmill is rapidly becoming an indispensable device for rehabilitation and conditioning for any age group.

Because of its unique safety, efficacy and comfort characteristics it is well positioned to benefit senior patients.

Approved uses include rehabilitation from any lower body disability, neurological gait training, geriatric conditioning, fall prevention and weight control programs.

The AlterG has been chosen by leading assisted living, skilled nursing facilities, hospitals and physical therapy clinics including; Life Care Centers of America, Steadman-Hawkins, UCSF, NovaCare, Stanford University, Harvard University and Walter Reed Army Medical Center.

Contact AlterG today to learn more:
T: (510) 270-5900 E: sales@alter-g.com

Or visit our website:
www.alter-g.com

Technical specifications

Weight reduction range:	0 - 80 % of BW
Weight reduction increments:	1%
Treadmill starting speed:	0.1 mph
Speed range:	0 - 12 mph
Speed increments:	0.1 mph
Incline range:	0 - 15%
Reverse speed range:	0 - 3 mph
User weight range:	90 - 400 lbs
AlterG shorts sizes:	XS, S, M, L, XL, 2XL, 3XL
Length x width x height	84" x 38" x 74"
Weight	750 lbs
Power requirements:	220V / 20A
Included warranty:	1 year
Available warranty:	Up to 5 years

Accessories

Video monitoring system for visual feedback and gait analysis. Thick rubber floor mat for floor protection and additional sound-proofing. Entry/Exit Step to facilitate access for limited mobility patients.

ALTER [®]
 Anti-Gravity Treadmill